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LASER BASED ANALYTICS

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Content

Vorwort	02-07
Preamble	
Einleitung	08-11
General introduction	
1 Laser-induzierte Fluoreszenzspektroskopie	12-26
1 Laser-induced fluorescence spectroscopy	
1.1 Einführung	14-17
1.1 Introduction	
1.2 Photolumineszenz	18-21
1.2 Photoluminescence	
1.3 Charakteristika der Fluoreszenzstrahlung	22-26
1.3 Characteristics of fluorescence radiation	
2 Hauptanwendungen der Wechselwirkung	28-55
von Materie mit Laserstrahlung	
2 Main applications of the interaction	
of matter with laser radiation	
2.1 Molekulare Lumineszenzspektroskopie und -analytik	30-35
2.1 Molecular luminescence spectroscopy and analytics	
2.2 MALDI-TOF	36-45
2.2 MALDI-TOF	
2.3 Opto-Akustik LIOAS	44-49
2.3 Opto-acoustics LIOAS	
2.4 LIPS	50-55
2.4 LIBS	

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LTB – 15 years of innovative laser techniques at Adlershof

You have all been linked, for many decades, with the scientific site of Adlershof. Why did you decide in 1990 to put an end to your scientific career and to establish your own company?

Dr. Matthias Scholz: We were all in our mid forties and had been working for 20 years at the institute for scientific equipment manufacturing. We felt, at this stage in our lives, that it seemed to be the ideal moment to put our many years of experience in the field of laser technology, into practice. The desire for independence was obviously of prime importance, the materialistic aspect was a secondary consideration.

One of the biggest hurdles in the founding of one's own company is the financial aspect. How was LTB Lasertechnik Berlin GmbH financed at first?

Claudia Scholz: As LTB Lasertechnik Berlin GmbH was founded according to GDR legislation, each partner had to supply 12,500 GDR marks for the company, which was, at the time, a very important sum. Three months later our nominal capital was devaluated from 2 to one in the course of the currency conversion. Overnight our 50,000 GDR marks became 25,000 German marks. According to the Federal Republic of Germany GmbH-legislation 25,000 DM were missing from the nominal capital. This extra capital had to be provided by the 4 founders and their families, with many sacrifices.

Behind every company there is normally a sustainable business model. Which business idea is LTB based upon?

Johannes Schurack: At the end of the eighties, laser technology was still a luxury which only well-known research facilities, universities and large affiliated groups could afford. Our business idea consisted of developing compact and economical laser systems so as to also allow institutes and companies with smaller budgets, to use laser technology.

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Founding a company in Germany, in view of all the administrative problems, requires an average of three months before every administrative barrier has been solved. How long did this period last for LTB and what was the most difficult part during the establishment of the company?

Dr. Stephan Mory: Our research works were the property of the Academy of Sciences within the former GDR. Each research domain was examined, regarding its relevance to the research environment of the Federal Republic of Germany, during the course of evaluation. Since this process was extended over many months, our original idea was, at first, uncertain. Fortunately, the appointed consulting company for the evaluation didn't consider our former research work as being marketable. Looking back on it now, it seemed a hazardous decision but it did lay the foundation of LTB. We could retrieve our patents and go ahead with the creation of the company. Another difficult point was the search for our own premises: at the beginning of the nineties, there were generally no commercial premises available in East Berlin. So we started looking and found a former police station that we were able to rent, after a lot of persuasion. Such a situation is not imaginable today, with the excess supply of commercial premises.

International companies such as Zeiss, Nikon, Canon and Bruker as well as renowned

research institutes are among your customers. With which strategy could a newcomer such as LTB, build up such a well known customer base?

Dr. Matthias Scholz: It is certain that our cooperation with the former Lambda Physik GmbH helped us. With a strong partner to back us up, doors opened which would normally have been closed for a small company such as LTB. Furthermore we have always tried to adapt to our customers' requests -

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even if we had to postpone our research ambitions from time to time. It is the customer's decision that determines whether it is going to be a success on the market, which doesn't necessarily mean high performing products equipped with all the latest technical finesse.

Technology orientated companies especially, often need years before making a profit for the first time. When did LTB make its first profit and how has the turnover developed since then?

Claudia Scholz: As we started with two marketable products and made a turnover from the beginning, we fortunately never had a deficit, but neither did we make a significant profit during the first ten years. We made our first profit in the year 2000, as we multiplied our turnover by tenfold, from 250,000 EUR to 2.5 million EUR, through lucrative industrial orders and increased it in the following years up to 6 million EUR. We are aware that these were exceptional years and that turnover development will not continue at the same growth rate. Since the investment policy of the industry has decreased noticeably, our turnover should level off to 4 million EUR in the coming years.

Taking a look at your product range shows that LTB is known for high performances: MNL- the most modern range of nitrogen laser, ELIAS- the world wide most efficient range of echelle spectrometer and ARYELLE- Butterfly- a new spectrometer generation, which combines two spectrometers into one. Where did your ideas for your product development come from?

Johannes Schurack: With our product range, it is similar to sport, only those standing on the podium are in the spotlight. Of course we also offer products which do not achieve top performances but which are, in quality, excellent. New product ideas arise from intensive dialogue with research facilities situated at Adlershof. Having the latest research trends directly on our door step, so to speak, we just need to pick them up and to convert them into marketable products.

According to Mr. Wolf, Berlin senator for economics, LTB is a prime example showing how, with a good idea, a lot of staying power and a significant engagement, a company can be set up successfully. Appraisal coming from a competent source; but surely there must also have been difficult phases?

Dr. Stephan Mory: As a small company from the east, it was particularly hard for us to

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maintain our hold on the market at the beginning. As an eastern-German company you had to be at least twice as good as the competitors. Entering the industrial markets was also challenging, as confirmed scientists, we were entering an entirely new domain and confronted with different requirements. In place of unit production a higher number of items were suddenly required. We wouldn't have been able to make such an effort, without Berlin public funding.

Within 15 years LTB progressed from being a regional company to a leading supplier for innovative laser systems. What were the milestones encountered?

Dr. Matthias Scholz: The most important milestones which paved the way to enter the industrial market at the end of the nineties were the improvements to our nitrogen lasers, to comply with international standards, and also the new development of our Echelle - spectrometers. With these two product ranges and with the manufacturing of the Aryelle series from the beginning of 2005, we progressed, so to speak, from the German Football League to the Champions' League.

But not wanting to rest on our laurels, a new generation of our nitrogen lasers will be presented at the laser exhibition in Munich, which should enable us to enter the field of mobile analytics.

More and more companies reduce jobs in order to improve their profit. You have provided 24 highly-qualified jobs. How does the present staff structure appear and how does an entrepreneur cope with social responsibilities within a tense economic situation?

Claudia Scholz: More than half of our employees are certified engineers, who come from the electrical engineering and physical and medical measuring technique fields as well as from informatics. Our female gender quote stands at 25 per cent, which is quite unusual in our sector. Furthermore we employ three trainees and three graduates. During the exceptional years of 2001 to 2003, we even had 35 employees. We always consider very seriously the social aspects when we have to make redundancies, because of the order situation, which is never easy for us. We always have the hope of getting the maximum number of people back on board when the economic situation improves.

Laser technology is among the fields which have the biggest growth potential. Which

Product training of Marubun in Japan by LTB, J. Schurack (1996)

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market position does LTB have within international competition and what are the strengths and chances of the company?

Johannes Schurack: LTB offers innovative products, efficient and economical laser systems for niche markets which are characterised by short product cycles.

We have acquired a leading position in a few fields, since, as a small company, we can react significantly more quickly to the current trend of the market. I am sure that we can gain further market shares in the field of laser based analysis methods, for example LIF or LIPS in the long term.

Being the site of the former science academy of the GDR, Adlershof has a long scientific tradition. What does the site offer today and what advantages are there for technology orientated companies?

Dr. Stephan Mory: Adlershof is a prime example of a successful synergy formation and efficient technology transfers. We have very close contact with numerous institutes. Intensive partnerships have often arisen from common R&D projects for the benefit of each partner. Adlershof is a successful model which should be used as an example for the whole of Germany.

Taking all these past experiences into account and looking back today, what would you do differently and what would you give as advice to company founders today?

Dr. Matthias Scholz: I wouldn't change anything, because it is mainly through mistakes that we have progressed. A decision made wrongly can be corrected at any moment; a decision which is not made at all, can have fatal consequences. As in football, key positions in a company need to be occupied by the right man or respectively the right woman. In order to be successful we need a centre forward who gives ideas, a strong midfielder for the practical implementation and a defender who pays attention to the fact that the ball doesn't go offside economically.

In past years the R&D investment of the German economy has continuously dropped. LTB regularly invests 10 per cent of the profits into new R&D projects. Do you manage better than other companies and why is R&D so important?

Claudia Scholz: The investment in research and development is important, in order to be able to react, with new products, to shorter innovation cycles. This aspect doesn't

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Hanover exhibition (1997)

only involve companies. Public investment can make a significant contribution and should therefore be enlarged. Because state aided R&D programmes are not only an investment for the future but are by all means, a profitable expenditure for the state. From the 2.2 million EUR subsidy that LTB has used, nearly 5 million EUR has gone back to the state, in the form of tax.

Technology orientated companies that want to expand their market position further in global competition, must always be ahead of their competitors with innovation. What are the present trends in your field and from which trends could a new LTB - product line emerge?

Johannes Schurack: In the analytical field, the trend is leading to economical solutions which must deliver an abundance of measuring data in a short time. In terms of rapidity we are well ahead with our ELIAS-series which enables 4000 measurements per second. Our nitrogen laser, with 30 per cent per series of measurements, is probably hard to undercut. In this context we are well prepared for the future.

Trade associations regularly criticise German bureaucracy and its complicated tax system. What can politics do to support medium size companies such as LTB?

Dr. Stephan Mory: The reform for company taxation is certainly going in the right direction. Lower tax rates for profits that are kept within the firm, will certainly have a positive impact on the investment policy of German firms. Furthermore, a new organization of the state loan programme would be desirable, up to now, R&D projects are only state aided up to the development of a prototype. However, there are very few state aided possibilities for the creation of marketable products.

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Introduction Interaction of matter with laser radiation –
induced and analyzed with lasers and measuring devices from LTB

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Specific signals

Interaction of matter with laser radiation in the wavelength range between near UV and near IR may in case of absorption result in the output of a specific signal from matter; apart from that matter leaves the interaction process unchanged. Examples of such specific signals are the laser-induced fluorescence (LIF), as well as the laser-induced pulsed temperature increase of matter and, resulting from that, a pressure wave, which is the basis of laser-induced optoacoustic spectroscopy (LIOAS). The change of refractive index, induced in the same thermic way, is the basis of time-resolved thermal lensing spectroscopy (TRTL).

Laser-induced fluorescence (LIF)

(or, more generally, laser-induced luminescence, which beside the normal fluorescence also includes delayed fluorescence and phosphorescence) may on the one hand give information on excitation and deactivation processes of the fluorescent species itself, moreover, the fluorescent species may act as probe in the investigation of complex systems, and last not least the fluorescence signal may be used to analyze classes of substances and to determine parameters of environment. To solve measuring problems from all three main application fields of LIF-spectroscopy and -analytics, LTB offers the complex measuring system LIMES. Based on LIMES, for special LIF-applications LTB offers especially adapted smaller measuring devices.

Laser-induced opto-acoustic spectroscopy (LIOAS)

may yield kinetic and calorimetric information about such short-living excited species, which deactivate more or less thermically, i.e. without emission. Therefore, this method is a useful counterpart to LIF-spectroscopy, especially also applicable to opaque or strongly scattering systems, e.g. complex biological systems. As laser source for this opto-acoustic spectroscopy, the wide range of pulsed lasers from LTB has proved successful.

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Light Scattering

Interaction of matter with laser radiation in the wavelength range between near UV and near IR may also result in light scattering. The scattered light may have the same wavelength as the exciting light (Rayleigh-scattering), or its wavelength may be shifted (Raman-scattering). Both types of scattering are used e.g. for optical distance measurement and remote ranging (LIDAR). Moreover, Raman scattering is the basis of a powerful spectroscopic technique. In both fields, the laser has brought substantial progress, and here especially the wavelength-tunable pulse laser radiation of LTB's nitrogen lasers/dye lasers is well suited.

Laser-LIDAR (light detection and ranging)

is based on measurement of the transit time of forward pulse and the back-scattered one; the method is useful not only for distance measurement and ranging of spatially fixed targets, but also of atmospheric clouds, clouds of smoke and exhaust, and generally for highly time-resolved distance measurements. In addition, especially with the variant of Raman-LIDAR, the back-scattered light can be detected spectrally resolved and this way the spatial distribution of species like CO₂, SO₂, N₂ can be measured. One special quality is the so-called resonance-LIDAR, it is based on laser-induced (resonance) fluorescence. In this case the backcoming signal may be much more intense than a normal scattered signal.

Laser- (resonance) Raman spectroscopy,

a method complementary to the infrared spectroscopy, is useful e.g. for the determination of inter- and intramolecular interactions and of molecular configurations, as well as for the energetic localization of optically “dark” excited electronic states. The LTB measuring system LIMES and its problem-adapted variants, respectively, are also suited for this spectroscopy.

Alterations of matter

Interaction of matter with laser radiation in the near UV, visible and near IR may at sufficiently high energy density also lead to irreversible alterations of matter. Under the gene-

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ric term “ablation” a lot of complex processes with widely varying aspects of application are classed. With respect to ablation the critical parameters of laser radiation are tunability of wavelength, possibility to get a sufficiently sharp focus, and problem-adapted pulsed excitation. Lasers from LTB with their broad variability of the said parameters offer tailor-made solutions for the various fields of application of ablation:

Laser treatment of materials

a) laser machining and laser-induced structural changes on the basis of local melting and/or vaporization, e.g. cutting, writing, drilling, welding, surface hardening, laser annealing. The wide-spread bow spans as far as e.g. the manipulation of biotissue, laser-microdissection, and laser surgery.
b) laser-induced chemical processing on pyrolytic or photolytic basis (thermochemistry/photochemistry), e.g. coating, doping and patterning of surfaces, synthesis of materials. The wide-spread bow extends to e.g. the laser based photodynamic tumor therapy (PDT).

Laser-induced breakdown (plasma) spectroscopy (LIBS)

Elementary analysis of solid materials by generation of microplasma with characteristic emission from excited atoms, ions or even molecular fragments. Possibility of increasing the LIBS sensitivity by combination with LIF. Special advantage of LIBS: Possibility of remote-site characterization. LIBS measuring systems from LTB: ARYELLE, ARYELLE-Butterfly.

„Matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF)“ mass spectrometry

Soft desorption and ionization of large biomolecules (e.g. proteins, peptides, oligonucleotides), which are deposited to suitable matrices to avoid fragmentation. Under the laser irradiation, these matrices vaporize, and the biomolecules are released to the gas phase with subsequent ionization and time-of-flight mass spectrometry. Also for MALDI-TOF, lasers from LTB have proved successful.

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These above mentioned characteristics of luminescence may reflect in an indirect, but sensitive manner also processes, which are in concurrence with the emission. Such processes either may belong to the intramolecular de-activation channels or can be intermolecular processes (Fig. 2). This way, the photoluminescence may deliver in a sensitive manner information on the emitting molecule itself (about its intra- and intermolecular dynamics) as well as information about its (micro)environment.


Fig. 2)
Manifoldness of processes, which may start from the first excited singlet state and may influence the fluorescence

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Thereof result the *main application fields* of molecular luminescence spectroscopy and analytics:

- Enlightenment of photophysics, photochemistry and structural relationships for the emitting molecule itself (in its particular environment);
- Use of the emitting molecule as a “fluorescent probe” in the investigation of complex biological, biochemical or physicochemical systems;
- Molecular emission as analytical method for the detection of multivarious classes of materials like e.g. drugs, environmental poisons, hydrocarbons, proteins, pigments; also the emitting molecule can be used as a molecular sensor for a definite physical/physicochemical parameter (e.g. pH).

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The availability of excitation by lasers has brought a real “revolution” in several of these fluorescence application fields in the access to information about materials and processes. The uniqueness and advantages, resp., of LIF extend in many dimensions; with respect to local dimensions it spans from the highest-resolution nanometer range to the kilometer range in remote sensing. Other important new fields of LIF application are based on laser-induced simultaneous or stepwise multiphoton absorption (Fig. 1).

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Photoluminescence of a molecule occurs as result of absorption of photons, by which an electron from a ground state orbital is transferred to an excited state orbital; the molecule is then in an electronically excited state.

Instead of energy uptake by photons, molecules can be transferred to an excited state also by other means of energy uptake and as a result, also show luminescence. Such possibilities are for example:

- excitation as result of chemical processes, which results in chemoluminescence,
- excitation via biochemical processes, which results in bioluminescence,
- excitation by frictional forces, which results in triboluminescence,
- excitation by ionizing radiation, which results in radioluminescence,
- excitation in an electric field, which results in electroluminescence.

Photoluminescence is the generic term for normal fluorescence, (usually simply called fluorescence), delayed fluorescence and phosphorescence.

Fluorescence is connected with the transition of an electron from a singlet excited state to the ground state without change in the spin orientation of that electron; this spin is in the excited state as well as in the ground state oriented antiparallel to a second electron in the ground state. Molecular energy states, in which all electrons belong to pairs with antiparallel spin are named singlet states. (In the figures, these states are labeled by symbol S). Therefore, fluorescence is a singlet →singlet transition². Such transitions are “allowed”, which means they take place with a high transition rate (typically 10^8 s^{-1}).

The initial state of phosphorescence is a triplet excited state (in the figures labeled by T), in this case the spin of the lifted electron is oriented parallel to a second (so-called unpaired) electron in a ground state orbital. Phosphorescence occurs as return of the lifted electron in that ground state orbital, connected with spin orientation reversal. Therefore, phosphorescence is a triplet→singlet transition. Such transitions are “forbidden”, they take place with small transition

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² In a more strict definition, fluorescence is a spin-allowed radiative electron transition; as such, it could occur also between two triplet states. In reality, emission from higher excited molecular triplet states is extremely scarce.

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³ ISC as well as phosphorescence are enabled by a so-called spin-orbit-coupling, which results in small mixture of triplet character to a singlet excited state and vice versa.

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
rate (typically $\rightarrow 10^4 \text{ s}^{-1}$). Direct triplet excitation from the ground state is, as a singlet→triplet transition, also forbidden; triplet states become populated on the “round-about route” of singlet excitation with a subsequent spin reversal process (intersystem crossing, abbreviated by ISC)³, e.g. $S_0 \rightarrow S_1 \rightsquigarrow T_1$ (see Fig. 3).

Usually, by absorption of photons, molecules are not transferred to pure electronic excited states, but simultaneously also vibrations become excited. These vibrations usually decay very rapidly, i.e. while the molecule is still in the electron excited state. In this case, photoluminescence starts from this vibrationless state. (In the figures, the labeling of the energy states (S;T) gives with the first index the number of the electron state, with the second index that of the vibrational state. Example: $S_{n,v}$ denotes the v-th vibronic state in the n-th singlet excited state).

With organic molecules, the starting level of direct fluorescence is usually the state $S_{1,0}$, independent of the initially excited state. There are only few exceptions from this so-called Kasha-rule; one of them is fluorescence from the second excited singlet in azulene. This rule is based on the rapid conversion of electronic excitation energy into vibrational energy (Internal conversion, IC), which may occur in higher excited singlet states (see Fig. 1, chap. 1.1). The same is true for higher triplets.

Fig. 3)
Schemes of excitation of phosphorescence and delayed fluorescence (two types) in comparison to normal fluorescence

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Beside the normal fluorescence also delayed fluorescence may occur. There are two different causes for this (see Fig. 3): Delayed fluorescence of “E-type” (named after its first observation with eosin) is based on the thermic excitation of a molecule from its vibration-less lowest triplet state to a higher vibronic state, from which ISC back to the singlet system is possible, with subsequent (delayed) fluorescence. Delayed fluorescence of “P-type” (at first observed with pyrene) may occur, if two identical molecules, each in the triplet state, interact in such a way, that one of them returns to the ground state and the other jumps into a higher excited singlet state. From there, it may relax in usual way by IC and then show (delayed) fluorescence. The said interaction process of the two molecules is called triplet-triplet annihilation (see Fig. 3).

Delayed fluorescence shows the same spectral distribution (color) like the normal fluorescence, but a different temporal behaviour: it delays more slowly.

Usually (by conventional spectroscopy) photoluminescence is excited by one-photon absorption. With lasers, excitation intensities are available, which render possible also *simultaneous* two- or multi-photon absorption. This means, photoluminescence can be excited with (instead of one photon of energy $E = h \cdot c \cdot \lambda^{-1}$) two photons of energy $E/2 = h \cdot c \cdot (2\lambda)^{-1}$ each, or with n photons of $E/n = h \cdot c \cdot (n\lambda)^{-1}$ (see Fig. 1). This way, the excitation wavelengths shift more to the red/near infrared spectral region. This may have advantages, e.g. for in-deep excitation in tissue.

Photoluminescence may be excited also by *stepwise* absorption of two or more photons, also this is a laser-based process. This way e.g. it is possible to detect ultraweak fluorescence from higher excited singlet states, what is impossible with direct one-photon excitation. Also this may have special advantages. Fig. 1, chap.1.1 shows the excitation scheme of such an anomalous fluorescence.

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Primarily, the fluorescence of organic molecules is characterized by its spectral distribution (color). To measure this fluorescence spectrum, the wavelength-dependence of the intensity distribution is measured at a suitably fixed excitation wavelength (e.g. at the maximum of the absorption band). With molecules in condensed phase, at room temperature usually this fluorescence spectrum consists of one band with a spectral width of several 10 nm. The reason for this remarkable width is that the emission line is exposed to spectral broadenings of different origin like rotations, impacts, and fluctuations of the molecule environment (e.g. of the solvent cage), even if the fluorescent electronic transition in question takes place solely between one vibrational state of a dominant vibrational mode in the terminal and final state each (e.g. $S_{1,0} \rightarrow S_{0,v}$ in Fig. 1, chap. 1.1).

In most cases, from $S_{1,0}$ there are further transitions with differing intensity ($S_{1,0} \rightarrow S_{0,\mu}$, $S_{1,0} \rightarrow S_{0,\nu}$). This so-called vibrational progression may be reflected by several sub-bands in the fluorescence spectrum, especially at low temperatures (Fig. 4, lower, left), or it can be smeared out (Fig. 4, lower, right). Also informative is the intensity distribution within such vibrational progression: If the intensity of the sub-bands is monotonously decreasing from the shorter to lon-

Dr. D. Leupold
LTB Lasertechnik Berlin


Fig. 4)
Potential energy diagrams (top) and related absorption/fluorescence spectra (bottom) of fluorophores without (left) and with (right) change of nuclear configuration as result of excitation. In addition, a different degree of spectral resolution is shown to illustrate the usual effect of low (left) and room (right) temperature.

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ger wavelenghts, this reflects almost identical nuclear configurations in the ground and excited state (Fig. 4, left). If on the other hand

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from shorter to longer wavelengths the intensity of the sub-bands goes through a maximum, this means that the relaxed nuclear configurations in the ground and excited state are different from each other, i.e. their potential surfaces are shifted from each other (Fig. 4, right).

In the simplest case, the fluorescence spectrum is the mirror image of the longest-wavelength (red-most) absorption band (see also Fig. 4). This mirror symmetry reflects, that the energetic differences (spacings) between the vibrational levels in the ground and excited state are more or less equal.

Usually, there is a spectral shift between the red-most absorption and fluorescence band, namely a red-shift of fluorescence. This so-called Stokes-shift is primarily caused by the FC-relaxation at differing electron density distributions as well as nuclear configurations in the ground and excited state. Due to influences from the environment, such differences in electron density distributions can result in even enhanced Stokes-shifts: for example, reorientation of a solvent cage as reaction of a changed electron density distribution after electron transition in the dissolved molecule results in an energetic stabilization of the excited state. Such solvent-caused red shifts of fluorescence (bathochromic solvatochromism) may be as large as several 100 nm.

The said effects like change of the solvent cage result in changes of the fluorescence properties only if they occur within the lifetime of the fluorescent excited state, e.g. if the solvent is sufficiently viscous.

In certain cases, also the excitation spectrum of the fluorescence may be informative. To measure it, the fluorescence intensity at a fixed wavelength (for example at the maximum of the fluorescence band) is registered while varying the wavelengths of a constant excitation intensity over the spectral range of the corresponding absorption spectrum. Differences between the latter and this excitation spectrum give e.g. for composed energy transfer chains direct information on the efficiency of intermolecular transfer.

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Measurement of the excitation spectrum of fluorescence also may be an important method to identify traces of organic impurities.

The fluorescence lifetime τ_F is the reciprocal of the sum of rate constants of all processes, which contribute to the deactivation of the fluorescent excited state; if these processes are exclusively intramolecular, then it holds

$$\tau_F = (k_T + k_{IC} + k_{ISC})^{-1},$$

where

k_T is the rate constant for radiative deactivation, which is determined exclusively by the electronic properties of the molecule.

k_T^{-1} is called natural lifetime or radiative lifetime, respectively. As a rule, it is the larger, the weaker the corresponding absorptive transition is.

k_{IC} is the rate constant of internal conversion (transformation of electronic in vibrational energy and radiationless deactivation of the latter)

k_{ISC} is the rate constant of singlet-triplet conversion ("intersystem crossing").

The more dominant k_T is in relation to both radiationless processes, the higher is the fluorescence quantum yield Φ_F ; it represents the ratio of the number of emitted and absorbed photons; it holds

$$\Phi_F = k_T(k_T + k_{IC} + k_{ISC})^{-1}$$

The definitions of Φ_{IC} and Φ_{ISC} are analogously; they represent the portion of excited molecules, which deactivate via internal conversion and via intersystem crossing, respectively.

All intermolecular processes, which become effective during the lifetime of the fluorescent excited state, contribute to the decay and diminish the fluorescence quantum yield, e.g. impacts with so-called fluorescence quenchers (called dynamical quenching); other types of such processes are indicated in Fig. 2, chap.1.2. On the other hand, there are also intermolecular processes, which influence the fluorescence quantum yield, but not the decay. These processes take place in another state than the fluorescent state, one example is the formation of quenching complexes in the ground state (static quenching).

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Another informative characteristic of fluorescence is their polarization: if molecules in an isotropic distribution are irradiated by monochromatic, linearly polarized light, then preferentially those molecules with transition dipole moment parallel to the direction of polarization become excited, which means that because of this so-called photoselection an anisotropic distribution of excited molecules occurs. If these molecules fluoresce, then also the fluorescence is anisotropic, i.e. is polarized preferentially in one direction. If I_p and I_s are the intensities of fluorescence parallel and perpendicular to the polarization of excitation, then the anisotropy A is calculated according to

$$A = (I_p - I_s) (I_p + 2 I_s)^{-1}$$

This anisotropy is maximal, if during the lifetime of the fluorescent state there is no rotational diffusion (so-called rigid systems) and no energy transfer.

Then it holds

$$A = (3 \cos^2 \delta - 1)/5,$$

with δ as the angle between the transition dipole moments of absorption and emission, respectively, which this way can be determined experimentally. For example, in case of parallel orientation of these moments it holds

$$A=2/5.$$

In case of perpendicular orientation (which may be the case for excitation to S_2 and fluorescence from S_1) it holds

$$A=-1/5.$$

In the other extreme case, if during the lifetime of the fluorescent state the excited molecules become completely isotropically distributed, it holds

$$A=0.$$

For most of the molecules, especially also for biomolecular systems, the fluorescence anisotropy is within these limiting cases; it decays during the fluorescence lifetime only partially. Here the static and especially the time-resolved measurement of the anisotropy (the latter realized by short-pulse excitation) gives informative access to rapid energy transfer processes as well as to movement, size and shape of macromolecules, measured e.g. via rigidly coupled fluorescent probes.

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Main applications
of the interaction of matter
with laser radiation

- 2.1 Molekulare Lumineszenzspektroskopie und -analytik
2.1 Molecular luminescence spectroscopy and analytics

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Molecular luminescence
spectroscopy and analytics

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The essential applications of molecular luminescence can be subdivided into the following three categories:

1) Luminescence spectroscopy is used to enlighten the photophysics, photochemistry (in the variety as shown in Fig. 2, cap. 1.1) and structural relationships for the emitting molecule itself (in dependence on the parameters of its micro-environment). In this respect, a comprehensive example is the (bacterio)chlorophyll fluorescence as a universal carrier of information about photosynthetic light harvesting, energy transfer, as well as energy dissipation processes in relation to structure-based excitonic couplings of pigments.

2) Use of the emitting molecules as fluorescent probes in the investigation of complex biological, biochemical or physico-chemical systems. Typical problems in this respect are e.g. the dynamics of polymer chains; lipid-protein interactions and light-induced changes in biological membranes; positions of bonds, conformation changes, denaturation, dynamics and geometries in proteins; interaction of cells with each other and with viruses. In these cases,

(i) the fluorescent probe may be itself part of the system under investigation (“intrinsic probes”), this is the (rare) ideal case, one example is tryptophan in proteins; or (ii) an extrinsic fluorescent probe may be covalently bound to the system under investigation, fluorescein derivatives with functional groups for covalent bonding to amino acids are an example, but also this case is rather rare because of the difficulty to synthesize macromolecules with such covalently bound, specifically fluorescing probes; therefore (iii) extrinsic fluorescent probes non-covalently associated to the system under investigation represent the dominating application. The site of solubilization of the probe within the system depends on the chemical nature of the probe as well as on the inducible specific interaction at that site. In the simplest case, the fluorescent probe is an indicator of a hydrophilic or hydrophobic environment.

In cases (ii) and (iii) a possible complication of interpretation may be caused by a local perturbation of the system under investigation by the introduced fluorescent probe; this per-

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turbation is determined by the shape (structure) and the size of the probe in relation to the probed region.

3) Use of molecular emission as analytical method. In this case, the target of analysis itself may be fluorescent (direct fluorimetric detection); examples are the detection of aromatic hydrocarbons, proof of proteins in blood serum and milk, detection of drugs, detection of chlorophyll. Most cases represent indirect fluorometric detection, e. g.

(i) reaction of the object of analysis with a species which results in the formation of a fluorescent compound (this is the standard case in liquid chromatography application in biochemistry and clinical chemistry, it is also the basis of several methods to detect ions);

(ii) fluorescence-quenching collision of the object of analysis with a fluorophor (useful for the detection of gases like SO_2 , H_2S , also for O_2 , dissolved e.g. in water or blood). A fluorophor, suitable for such a selective detection of a species, is called molecular fluorescence sensor. Often the complete measuring device for selective detection based on such a sensor is also called fluorescence sensor, e.g. for O_2 .

(iii) There are also such molecular fluorescence sensors for certain system parameters, like pH. (This special case, represented e.g. by coumarines, is often called directly “fluorescent pH indicator”).

For these applications of molecular fluorescence spectroscopy especially in biological systems, a series of standardized techniques has been developed, e.g. FRET for the determination of intermolecular distances, based on the fluorescence measurement as result of Förster Resonance Energy Transfer, FLIM, the Fluorescence Lifetime IMaging of lifetime-sensitive fluorescent species, FRAP for the determination of lateral diffusion coefficients, e.g. of lipids or proteins, based on the laser bleaching of molecular probes (Fluorescence Recovery After Photobleaching).

Many of these applications of molecular fluorescence spectroscopy and analytics either has become possible only since the availability of lasers as excitation source of fluorescence, or lasers essentially improved their content of information. Some of these highlights are time resolution in the femtosecond range,

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remote sensing in the kilometer range, local resolution in the nanometer range as well as non-invasive optical biopsies. (To this category belongs also the subsequently described diagnostic application of a tissue fluorescence excited via a femtosecond laser-based stepwise two-photon absorption).

Application example

Femtosecond laser-based analysis of melanin fluorescence (Early diagnosis of black skin cancer)

Up to today, fluorescence-based methods do not play any role in the differential diagnosis of malignant melanoma. Now, a femtosecond laser-based equipment is ready for clinical testing, which will allow fluorescence spectroscopic identification of even smallest areas of tumorous skin tissue. The method is based on a novel mechanism of stepwise two-photon absorption of melanin.

The role of melanin in human skin may be characterized as a two-edged sword: On the one hand, melanin is responsible for pigmentation of skin and its parasol-like protection. But, on the other hand, melanin also plays a key role in the malignant degeneration of healthy skin tissue. Therefore it seemed well-grounded to look for a melanin-based diagnosis method for malignant melanoma. Until quite recently however, such investigations were not successful, especially with fluorescence.

There are two reasons for it: firstly, the fluorescence of melanin is extremely weak; only about each tenthousandth of the absorbed photons is re-emitted as fluorescence, i.e. the fluorescence quantum yield is about 10^{-4} , the predominant portion of the absorbed excitation energy is transformed to thermal energy. Secondly, the fluorescence spectra of healthy tissue and of malignant melanoma (when measured in the conventional way by one-photon excitation, Fig. 6a) differ only very slightly (Fig. 5) and moreover, the fluorescence of other endogeneous fluorophors with considerably higher fluorescence quantum yield like NADH, tryptophane, flavins may additionally mask this



Fig. 5) Fluorescence spectra of excised samples of human skin tissue, obtained by conventional one-photon excitation at 400 nm (Fig. 6a). Blue: normal healthy tissue; red: malignant melanoma (nodular melanoma, Clark level 4, tumor thickness 1.7 mm according to Breslow).


Schematic representation of
fluorescence excitation

Fig. 6a)
in the conventional manner by
absorption of one photon


Fig. 6b)
by stepwise (subsequent) absorption
of two photons of half the energy
(twice of the wavelength) as com-
pared to case a).

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difference. Experiments with markers did not improve the situation.

In summary, till now the general criteria of a clinical method for early detection of malignant melanoma could not be fulfilled.

- Such criteria, adapted to the aim of melanoma detection by fluorescence, are
- selectivity of the method (especially also in-vivo, this means exclusion of any disturbance by endogeneous or exogeneous fluorophores), and
 - the malignant degeneracy parameter must change the selected fluorescence parameter specifically and significantly.

The solution of the problem, i.e. a new approach to the melanin fluorescence in-vivo, could be found only recently by making use of the unusual absorption behaviour of melanin. Its spectrum between the UV and the near IR shows a monotonously decreasing shape, whereas other organic compounds typically show absorption bands in distinct spectral regions. Due to this special absorption, it is possible to excite melanin in a spectral range, where no other components of skin tissue absorb and therefore do not fluoresce. Such spectral range e.g. is at and around 800 nm. The energy of one 800 nm-photon is too small to induce the melanin fluorescence. But thanks to an intermediate energy level melanin can absorb a second 800 nm-photon, and via this stepwise two-photon absorption melanin reaches the fluorescent level (Fig. 6b).

The fluorescence spectrum of the this way selectively excited melanin reflects specifically and significantly a malignant degeneration of skin tissue. The significant change of the spectral shape of the fluorescence is shown exemplarily in Fig. 7: the malignant degeneration shifts the maximum towards the red. Investigation of a broad range of other skin diseases like basal cell carcinoma, fibroma, keloid, balanitis, squamous cell carcinoma showed, that this spectral red shift of the fluorescence occurs exclusively with malignant melanoma. The cause of this shift is a changed ratio of pheomelanin to eumelanin: In the process of malignant degeneration this ratio changes in favour of pheomelanin. This

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could be concluded from laser-fluorescence spectroscopic investigations of prepared different melanin species.


Fig. 7) Fluorescence spectra of excised samples of human skin tissue (the same as in Fig. 5), obtained by stepwise two-photon excitation by an 800 nm-photon emitting femtosecond laser (Fig. 6b). Blue: normal healthy tissue; red: malignant melanoma (nodular melanoma, Clark level 4, tumor thickness 1.7 mm according to Breslow).

The investigations also showed, that in the process of malignant degeneration a second fluorescence property changes: The decay time of the fluorescence of the malignant melanoma is distinctly shortened as compared to the fluorescence decay from healthy tissue. Also this is caused by the changed pheomelanin/eumelanin ratio. The fluorescence lifetime of pheomelanin is shorter than that of eumelanin.

With the demonstrated significance as well as specificity both the suppositions of a reliable early diagnosis of malignant melanoma by fluorescence diagnostics are fulfilled.

As light source for the stepwise two-photon excitation of the fluorescence of skin tissue a femtosecond-pulse laser is used. It delivers the necessary high intensity during the short pulse duration, which in one realizes the optical preparation of a sufficient concentration of fluorescent species and avoids tissue damage. Moreover, due to the ability to focus the laser radiation early diagnosis on very small skin areas is possible.

After successfully finishing the basic research with laboratory equipments at stationary femtosecond lasers, a first mobile, intensified femtosecond fiber laser-based device for fluorescence diagnostics of malignant melanoma (MLS 200, LTB Lasertechnik Berlin) is tested for clinical use.

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
Introduction

Mass spectrometry is an analysis technique to determine by means of a mass spectrometer the molecule mass of free ions in the high vacuum.

A mass spectrometer (Fig. 8) consists of an ion source where a beam of gaseous ions is formed from a sample, a mass analyzer that separates the ions regarding their mass-charge-ratio (m/z), and a detector to deliver a mass spectrum. The latter shows which ions in which relative quantities are formed.

Since the 1970's laser devices were employed in mass spectrometric analysis to achieve a direct desorption of intact molecule ions from condensed gas phases by a suitable primary excitation. During these initial attempts a thin layer of the samples was placed on a metal surface to be irradiated by a pulsed laser. However, mass spectra resulting from these approaches were of less intensity, as the laser energy was sufficient enough to cause molecular fragmentation of the sample molecules. Therefore the laser desorption ionization mass spectrometry (LDI-MS) was of less practical use for the analysis of bio molecules since only ions of a molecular weight below 1000 Da could have been detected.

This situation changed entirely as M. Karas and F. Hillenkamp from the University of Münster inspected the correlation process between UV laser irradiation and organic molecules. They figured out that embedding the sample into suitable matrix material, consisting of small organic molecules, the laser light was strongly absorbed and

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thus not only higher intensities of the analyte were obtained but also sample fragmentation was reduced to a minimum. This method denoted MALDI-MS (Matrix Assisted Laser Desorption/Ionization) allows now detecting large molecules such as proteins.

Analyzers associated with MALDI experiments are Time-Of-Flight instruments where the mass determination in the high vacuum area is performed by a very precise measurement of the period of time after the acceleration process of the ions in the source and their impact on the detector, this combination is called MALDI-TOF.

Basics

The first step of mass spectrometry is ion formation. Ionization of the analyte molecule is performed by reception or loss of an electron. This process can be initiated by electron impact (EI), bombarding the sample with atoms, or ions (Fast Bombardment Ionization, FAB), or by photons (Laser desorption/ionization, LDI). A further ion formation process is denoted Electro Spray Ionization (ESI) where solved sample is sprayed into an electrical field.

Ion formation in the MALDI process works as follows: mixing an analyte on a metal sample probe with a suitable matrix compound of a 1.000 up to 10.000 times molar excess absorbing at the used laser wavelength, causes a cocrystallization process of both matrix and analyte material after evaporation of the solvent. The incorporation of the sample molecules into the lattice structure of the matrix is supposed to be the precondition that the laser desorption/ionization process will work.

The crystallized surface of the prepared sample is then exposed to an intensive, short wave laser pulse in the high vacuum area inside the ion source of the mass spectrometer. The coupling of the energy, which is necessary for the ions, is performed on UV radiation by resonant excitation of the matrix molecules, e.g. into the π -electron system of aromatic compounds. Theoretical

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calculations let suppose that the excitation energy being stored in the matrix molecules relaxes in extremely short time periods into the solid-state lattices causing here a strong distortion of the expansion. This process is followed by transition into a phase in which a part of the solid-state surface is vaporized explosively long before a thermal balance occurs. In this step matrix molecules and sample molecules as well are released into the gas phase (Fig. 9). Obviously the internal degree of freedom of the molecules sharing this process is low enough that even thermal labile macromolecules such as proteins endure this process. However this is only valid as long as a limited range of irradiation power on the sample between 10^6 and 10^7 W/cm² is concerned. Too a high radiation destroys the sample to a great extent. A lot of experiments indicate that the matrix is also considered an agent ionizing the sample molecules. Accordingly photo ionized, radical matrix molecules cause by transferring protons a high yield of electrically charged sample molecules. An electrode, which is mounted some mm apart opposite to the sample position, is used to generate an electrostatic field in the range of some to some 1000 V/cm. Depending on the polarity positively or negatively charged ions are accelerated from the sample surface towards the analyzer.


Fig. 9)
Principle of the MALDI-Process

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Laser devices employed in MALDI experiments are UV lasers, such as Nd-YAG lasers (Yttrium-aluminum-garnet-crystals doped with neodymium) at pulse durations of about 5-15 ns and at wavelengths of 335 nm (tripled frequency) or 266 nm (quadrupled frequency) respectively or N₂-lasers operating at a wavelength of 337 nm with a pulse duration of 1-5 ns. Concerning the infrared range Er:YAG-lasers (Yttrium-aluminum garnet-crystals doped with erbium) are at the disposal working at a wavelength of

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

$2,94 \mu\text{m}$ with a pulse duration of 90 ns. The beam is focused by a suitable ion optics onto the sample inside the ion source to a diameter of $\leq 150 \mu\text{m}$.

A special point of MALDI is the matrix preparation. The principle of MALDI assumes the insulation of analyte molecules in solid crystallized matrix material. As mentioned before one of the requirements to the matrix molecules is absorbing the irradiating laser photons. Additionally the analyte molecules must be incorporated into the lattice structure of the matrix in a manner to neutralize their intermolecular correlation. However the property of co-crystallization is associated only with a few number of compounds absorbing at the UV range. For example, from all the possible dihydroxyderivates of the benzoic acid only the 2,5 dihydroxybenzoic acid can be used.

The second step of mass spectrometry is ion separation. Using the Time-Of-Flight method, an electro-static field accelerates ions formed during a short laser pulse inside the source to a kinetic energy of some keV. After leaving the source the ions pass a field free drift region in which they are separated due

*Fig. 10) and Fig. 11)
Reflector flight tube*

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to their m/z ratio. This takes place because at a fixed kinetic energy ions with different m/z values are accelerated in the ion source to different velocities. Knowing the acceleration voltage and the length of the drift region the m/z ratio can be determined by measuring the flight time.

Fig. 10 and 11 illustrate the principle how a TOF mass spectrometer works in linear and

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reflector mode. Ions formed by the laser pulse of the same charge but with different m/z values are accelerated to different velocities. Large ions with high m/z values strike the detector at a later moment than small ions. Ions of the same mass start with a certain energy spread. This affects the peak width to a certain amount. Applying a reflector the influence of the initial energy on the flight time can be compensated at the location of the reflector detector.

Associated with a TOF instrument the relation of molecule mass and charge number (m/z) is proportional to the square of the flight time t . Thus the related mass can be determined from the measured flight time. Calibration is performed with well-known reference masses. The range of typical MALDI flight times is between few μs and some 100 μs . The drift regions are typically 1 - 4 m long.

An important process which may happen during the acceleration process, or when passing the drift region is a decay of molecule ions. The process of metastable decomposition in the area of the drift region was called PSD (Post Source Decay).

Metastable decomposition of the analyte ions happens by absorbing high energies during MALDI by collision excitation processes inside the matrix cloud immediately after the laser pulse has hit the sample. The efficiency of the metastable decomposition depends on the density of the matrix cloud after the laser shot and the height of the kinetic energy making analyte ions and matrix molecules collide. The density of the matrix cloud depends on the laser intensity and the temperature the matrix sublimates. The applied acceleration voltage defines the collision energy. Normally during metastable decomposition a charged fragment ion is formed and a neutral fragment molecule as well.

As a rule the analysis of PSD ions delivers useful information concerning the structure of the original molecule ion. For example, PSD analyses related to proteins benefit in sequencing techniques of peptides, since decompositions mainly occur along the peptide chain.

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The third step of mass spectrometry is ion detection. Generally, ion detection systems consist of a Secondary Electron Multiplier (SEM). The assembly works in a manner that impacting particles such as ions and neutrals produce secondary electrons, which are available as an electrical signal at the output of this assembly for further processing.

Modern detectors used in TOF mass spectrometric analysis are designed as micro channel plate detectors. Such a device is a


Fig. 12)

Scheme of a typical instrument arrangement of a MALDI measurement

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porous solid core assembly with millions of tiny holes ($\sim 5\text{-}10\ \mu\text{m}$, length $0.5\text{-}0.8\ \text{mm}$), which are coated inside with a semiconductive layer. The holes are denoted micro channels. Everyone of them works as an electron multiplier, independent of the adjacent ones. To obtain a maximum electron yield all these channels are electrically connected in parallel. Micro channel plate detectors benefit in detecting small to large input signals without changing the multiplier characteristics, thus keeping saturation effects to a minimum. The multiplier output voltage is guided to a transient recorder to be digitized and delivered to a PC (Fig. 12).

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Application example

The applications for the MALDI TOF/TOF instruments are in the life science business. In the today's major application, the protein identification, the Ultraflex II (Fig. 13) TOF/TOF is used to determine the weight of proteins out of complex mixtures. This very accurate mass determination which is equivalent to 1 gram out of a ton allows identification of big molecules. This is achieved by searching huge public protein databases and identifying the compounds just by knowing their mass. The goal behind this business called proteomics is to understand the cell biology on a molecular level, to finally discover pathways and develop corresponding drugs with a minimum of adverse reactions. The detection and use of biomarkers to detect diseases in an unprecedented early stage, just likewise important, is done by this approach in a field called "clinical proteomics".

Fig. 13)
MALDI-TOF instrument
„Ultraflex II“

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Introduction

In the industrial process control, an increasing importance is given to the non-destructive layers and surface testing. The life time and reliability of many modern components are essentially determined by the controlled characteristics of the coating. For the optimizing and examining of the layer characteristics, reliable and informative layer measuring methods are needed.

A series of imported process like micro and nano-hardness tester, scratch tester, impact and cavitation tester are available. The applicability is limited by the layers thickness, hardness, adhesiveness, surfaces roughness or porosity.

The laser acoustic test method LAwave from the " Fraunhofer - IWS" in Dresden offers a good complement to these testing method and opens further on far reaching new possibilities.

Wideband surface acoustic waves are excited in materials via sub-nanosecond laser pulses. Wave dispersion reacts very sensitively to coating and surface modifications.

The laser acoustic test method LAwave measures the dispersion velocity of the surface acoustic wave depending on frequency and can thereof determine material characteristics like modulus of elasticity of nanometer layers or the depth of damage layers within semi-conductor materials. The layers can be hard like diamond and soft like Polymers. It can also be measured on rough surfaces.

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
Measuring method

Surface acoustic waves are elastic vibrations which are propagated along the material surface. Their amplitude is at the maximum height at the surface and decays exponentially in the inner material. Their penetration depth is equivalent approx. to the wavelength. If the frequency increases, the penetration depth of the wave reduces. Three different cases showing the dispersion velocity of surface acoustic waves depending on frequency are described in Fig. 1: uncoated silicon, silicon with a diamond-like carbon film which is stiffer than substrate material, and silicon with a polymeric film which is considerably softer than substrate. In the case of the uncoated silicon, the velocity of acoustic propagation is constant because waves not depending on frequency and therewith not depending on penetration

Fig. 14)
Phase velocity c of surface acoustic waves on coated and uncoated silicon depending on frequency

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depth propagate themselves in the same material. This changes when the base material is coated. The diamond-like carbon film generates a rise in the curve as with higher frequency the wave penetrates not as deep in the material and therewith is more affected by the stiffer layer than by the substrate material. The layer accelerates the wave dispersion with increasing frequency.

The polymer layer generates the opposite effect. The wave velocity sinks with the frequency, as the coating material is considerably softer than the substrate. If the curves are extrapolated to zero point of the frequency they meet each other in the same point, the acoustic velocity of the substrate. This has to be put down to the fact that waves with decreasing frequency penetrate deeper and deeper in the material so that the layer exercises no more influence. Fig. 14 suggests that the ascend of the curve contains information on

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
Fig. 15)
Schematic representation of the
laser-acoustic equipment.

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the layer characteristics, on its modulus of elasticity, on its density and on the thickness of the film. From the position of the zero point, the modulus of elasticity of the substrate can be determined. The curves in Fig. 14 are called dispersion curves. They are measured with the laser acoustic measure system LAwave (Fig. 15) and their material characteristics are thereof calculated.

Laser pulses of a nitrogen pulse laser (pulse duration: 500 ps, pulse energy: 400 μ J, wavelength: 337nm) excite the acoustic surface waves which are detected by a piezoelectric transducer. The phase spectra are calculated via Fourier-transformation from the laser acoustic signals, taken in various distances of the detector from the laser focus line, from which the phase velocity C depending on the frequency f will be determined. The material parameters are deduced by fitting a theory to the measured curve. Such a dispersion curve can also be calculated with a theory based on the solution of the wave equation for the boundary conditions from tension and shifting at the surface and at the bordure area for layer and substrate. The theoretical curve is fitted to the measured curve via non linear regression calculation to determine material parameters.

The test needs only a few minutes. It can be conducted on plain or cylindrical surfaces. There are no further requirements to the specimen geometry. This measuring method is the more sensitive towards thin layers the lower the acoustic attenuation in the material is. Ultra thin layers can be examined on mono-crystal and fine-grained material like silicon or coated steel. The layers should not be much thinner than 100 nm for coarse-grained materials like stainless steel and aluminium.

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Application examples

Testing Ultra-thin Diamond-like Carbon Films

Evaluating ultra-thin films of a few nanometers thickness is still a challenge for mechanical testing. The laser-acoustic technique was demonstrated to have the capability of measuring the Young's modulus of diamond-like carbon (DLC) films with thickness of 3 nm deposited on silicon wafers.

Fig. 16 shows an example for a laser-acoustic results obtained for two DLC-films with thickness of 3 and 16 nm, respectively. The films were deposited by the pulsed arc technique on silicon. The measurement provides the disper-


Fig. 16)
Laser-acoustic measuring result for
ultra-thin diamond-like carbon
films deposited on silicon

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die Dispersion curve as phase velocity depending on frequency. The effect of the 3 nm DLC film causes a distinct rise of the dispersion curve. The Young's modulus is deduced by fitting the theoretical curve. The fit yielded a Young's modulus of 697 GPa. To verify the effect of the measuring uncertainty, the measurement was repeated 10 times. The film modulus scattered within the limits of " Δ " $E \pm 40$ GPa.

For diamond-like carbon, Young's modulus was shown to correlate with the sp^3 -fraction of the atomic bonds that also determines the hardness of this film material. That is the reason why Young's modulus correlates with hardness ($H > E/10$). The high film modulus proved the high quality of the ultra-thin diamond-like carbon film deposited by a new plasma source consisting of a high current pulsed vacuum arc evaporator combined with a sectioned filter unit]. The dispersion curve of the 16 nm film descend revealing the low elastic modulus of this film compared to the 3 nm film. The fit yields a film modulus of only 107 GPa. The film was deposited at higher substrate temperature favoring the graphitization of the carbon

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film. The tests prove the laser-acoustic test method to be able to indicate clearly the different mechanical behavior of these films and to provide reliable material parameters for the film material.

Testing Damage Layers in GaAs Wafers

Slicing GaAs wafers from ingots produces damage layers that have to be completely removed by polishing. A non-destructive control of the machining process is desirable. The laser induced surface acoustic waves has been demonstrated to be a method for characterizing the finishing state of such surfaces.

Fig. 17 presents the laser-acoustic results for a wafer stepwise polished after having been sawn from the GaAs ingot. For the as-sawn wafer (polishing depth $d = 0 \mu\text{m}$), the velocity drastically decreases with increasing frequency.


Fig. 17)
Laser-acoustic results for GaAs wafers with different polishing depth d

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This suggests that a layer with a remarkable lower velocity covers the surface. High density of micro-cracks must be assumed causing such a considerable reduction of the stiffness of the surface layer. Polishing the samples results in more gradually descending dispersion curves. The velocity is nearly constant for a polishing depth of $10 \mu\text{m}$. This means that all waves independent on the velocity and therefore consequently independent on their penetration depth propagate through the same material.

The results enable to use the slope for determinations of damages depth.

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Introduction

Laser-induced plasma spectroscopy, LIPS, is a method utilizing laser ablation and the subsequent atomic emission from the plasma for elemental analysis. Besides the acronym LIPS, today other designations such as LIBS (laser induced breakdown spectroscopy), LAOES (laser ablation optical emission spectroscopy), or LSS (laser spark spectroscopy) can be found in the literature.

Laser ablation is at present the only analytical method that offers direct sampling from any kind of material without sample preparation. So LIPS allows a multielement analysis of virtually all types of materials (gas, solids, liquids) through atomic emission spectroscopy. Today's availability of reliable and less costly laser sources and improved spectrometers and detectors permits a rapid, on-line, and in-situ analysis with LIPS. This makes LIPS especially attractive for all kind of process analysis and environmental screening and monitoring.

*Prof. Dr. U. Panne
BAM
Federal Institute for Materials
Research and Testing*

Basics

Fig. 18 displays the typical temporal sequence of a LIPS analysis. The laser excites vibronic and electronic states, the energy of which is converted to kinetic energy and finally results in the formation of a weakly ionized expanding plasma. The observed modification of the surface is beyond a monolayer and non-linearly related to the number of excited states. Further, an irra-


Fig. 18) Temporal sequence of a LIPS analysis.

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dianze threshold is observed for laser ablation. Due to the dominating unspecific emission from bremsstrahlung and recombination, the early phase of plasma emission cannot be used for atomic spectroscopy. So after a certain delay, relative to the plasma ignition, the characteristic atomic or ionic emission lines can be observed. The line intensity is integrated for a defined time and then utilized for analytical purposes. The observed emission is directly related to the concentration of the analyte under study in the sample resp. in the plasma. Typically, electron densities between 10^{16} - 10^{18} cm^{-3} and temperatures of 6000-23000 K are observed after 1-5 μs delay at 1 atm of air.

A distinct advantage of LIPS is not only the conceptual simplicity, but also the versatility and simplicity of the experimental set-up. Most present LIPS systems are based on


Fig. 19)
Echelle spectrum from a tungsten mineral (inset: 128th order, 302.0-304.8 nm); adapted from ref. (Florek et al. 2001)

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either a Q-switched Nd:YAG-laser or a gas laser such as excimer and nitrogen lasers for plasma ignition. Both the laser for plasma ignition and the plasma emission itself can be delivered via fiber optics, which allows a remote sensing over several meters. Another option is a direct sensing via telescopic observation and plasma generation. Conventional monochromators or spectrographs can yield either a wide spectral coverage or a high spectral resolution ($\lambda/\Delta\lambda > 10000$), which results in limited applicability for multielement analysis in complex matrices. For simultaneous detection of multiple elemental lines over a large spectral range, echelle spectrographs in combination with CCD detectors are the systems of choice. Echelle systems offer the necessary spectral coverage and spectral resolution for minimizing spectral interferences in complex

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matrices. Fig. 19 demonstrates the wealth of analytical information available from an echelle system, when complex matrices are studied.

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Application examples

Solid sample analysis is probably the most successful area of application for LIPS. Due to the fact that metals have been the preferential samples for all kind of fundamental and analytical studies with LIPS since the first investigations of laser ablation, alloys and steel samples have been investigated in detail either in laboratory-based studies or


Fig. 20a)
Experimental set-up for remote analysis of glass melts via LIPS (adapted from (Panne et al. 1998))

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process analysis work. Usually, the LIPS detection limits rival the classical spark analysers and additionally permit the analysis of non-conducting samples as slags.

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Another very popular area of LIPS analysis is the fast screening of natural materials such as soils, sands, and waste materials (e.g. sewage sludge). Although several problems were reported due to the heterogeneity of the matrix both spatially and chemically, usually detection limits in the order of 10-100 ppm are obtained if the sample can be conditioned through drying, mixing with a reference material, or pressing of pellets. Other materials of industrial importance

Fig. 20b)
Remoste sensing with LIPS through ignition of a laser plasma on a glass melt at 1200°C

Fig. 21)
LIBS analysis of recycled thermoplasts from consumer electronic products.

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investigated by LIPS are glass and ceramics, which can be also analysed in remote fashion. Fig. 20a shows the experimental set-up for a remote LIPS analysis of hot glass melts, while Fig. 20b demonstrates the plasma formation on a hot glass melt at 1200°C. Further LIPS applications can be found for ores and coals, paints and recycling materials such as plastics with different kind of fillers and pigments. Fig. 21 reveals a LIPS set-up installed at a polymer extruder to monitor the extrusion of recycled thermoplasts from consumer electronic products for toxic heavy metals.

Promising appears also the microanalysis of heterogeneous materials. Through special beam homogenization, a laser can be focused down to a spot size of 3-10 µm. LIPS provides not only spatial resolution, but offers also some depth resolution through multiple ablation on a single spot. In comparison to other beam techniques, LIPS needs only minor sample preparation and no vacuum. This approach was demonstrated for different heterogeneous and layered materials such as paper coatings, catalysts, or microinclusions in steel.

The detection limits for the direct analysis of dissolved analytes in water or colloidal materials in water significantly decrease (in the order of 0.1%) due to the fast deexcitation of the excited species. Not surprisingly, the most successful approach for analysis of liquids via LIPS is the ablation from the surface. In this case the plasma is expanding into a buffer gas, which permits ppm detection limits.

Analysis of gases and aerosols with LIPS was mainly studied for volatile metals especially mercury and halogenated compounds with detection limits in the ppm range. For the analysis of aerosols, two principal experimental set-ups, can be realized: (i) a direct, on-line detection in the focal volume or (ii) sampling of the aerosol on filter or impactor. While the first approach permits an on-line detection, it is hampered by the low sampling volume and the corresponding decreased sensitivity (detection limits in the µg m⁻³ range). Due to the threshold dependence of the plasma ignition, the pulse energy can be adjusted in such a way that only with an aerosol in the focal volume a plasma is generated. This allows a dedicated single pulse analysis which can improve the detection limits considerable. The analysis

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on filter samples allows to integrate larger volumes of air so that analysis of lower aerosols concentration (in the ng m^{-3} range), e.g. from remote areas, is possible. Aerosols were analyzed via LIPS not only for industrial hygiene and anthropogenic emissions, but also for climate change related studies

In summary LIPS represents for some of today's analytical problems the only reasonable solution to in-situ and/or elemental remote sensing. In future, the number of specialized applications will probably increase in the future due to the minimum sample preparation, the reduced cost for a single measurement, and the potential for an extensive automation. The application of solid sample analysis, especially microanalysis techniques for nanotechnology or ultra-small volumes, will grow due to inherent blank limits in ultratrace analysis.

Due to the fact that LIPS is technology-driven, advances in the laser technology and optoelectronics will have direct impact on the design and miniaturization of LIPS systems. Applications in process analysis will further profit from fiber-optical systems delivering the laser to the point of demand and guiding the emission to a centralized detection system. The combination with other types of spectroscopy providing molecular information such as Raman spectroscopy will open further applications for LIPS.