

**LIGHT. PRECISION. ANALYTICS.**


# LIBSpector

*Benchtop laser class 1 sample chamber for spatially resolved LIBS analysis*



## Key features

- For solid, liquid and gaseous substances
- Qualitative and quantitative analysis
- Motorized XYZ stage with joystick
- Pilot laser and video monitoring
- Built-in LIBS laser head optional




LIBS spectrum of steel measured with LIBSpector and ARYELLE 200 spectrometer.

The LIBSpector is a compact and intuitive to operate sample chamber for the qualitative and quantitative analysis of samples by means of laser-induced plasma spectroscopy (LIBS). In combination with a high-resolution echelle spectrometer from the ARYELLE and DEMON series and a suitable laser for plasma ignition the LIBSpector is part of our modular LIBS analyzer called LIBSlab – suitable for universal application in science as well as in industry-related quality control.

The beam of the LIBS laser, whose head can be built-in the chamber housing, is directed to the sample via telescope optics and generates a light emitting plasma. The plasma light is guided to the spectrometer via mirror and fiber optics. Sample mapping is provided by an integrated motorized XYZ stage that can easily be controlled either with the supplied XYZ joystick or by software.

All measurements can be monitored in real time with the installed high-resolution CMOS observation camera. The focus of the plasma-generating laser is at the intersection point of the pilot laser beams. Precise sample positioning and focusing is required to ensure accurate LIBS measurement results by simply adjusting the XYZ stage. The optimized optical setup of the model  $\mu$ -LIBSpector provides a highly reduced laser spot size on the sample surface. Several sample holders for solid, liquid and gaseous substances provide universal application capability and can be adapted to individual customer requirements.

## LIBSlab = modular system



The LIBSpector comes with a laser class 1 housing and is equipped with safety interlock and laser protection window for safe observation. No additional laser safety precautions are therefore required at installation site. Sample-specific dust and particles resulting from laser ablation are extracted by an externally connected exhaust unit.

## Applications

Benchtop laboratory equipment for scientific and industrial applications

With the LIBSpector you get a technically extensively equipped and functional LIBS sample chamber in modern benchtop design. It provides excellent conditions for a flexible and safe use of LIBS technology for qualitative and quantitative analysis in the scientific and industrial sectors.

## Specifications LIBSpector series

	<b>LIBSpector</b>		<b>μ-LIBSpector</b>
<b>Features</b>	Standard benchtop device for universal laboratory applications in science and industry		Optimized optical setup providing a very small laser spot on the sample surface
<b>LIBS Laser * (typ.)</b>	Wavelength	nm	1,064
	Pulse energy	mJ	50
	Frequency	Hz	20
	Laser type		Nd:YAG
	Laser spot size on sample	μm	∅ = 300
	Built-in laser head		optional
			∅ = 50
			optional
<b>XYZ stage</b>	Travel range	mm	X = 160, Y = 100, Z = 50
	Resolution **	μm	5
	Repeatability **	μm	± 10
	Velocity (typ.)	mm/s	v <sub>x</sub> = 25, v <sub>y</sub> = 25, v <sub>z</sub> = 5
<b>Components</b>	Observation camera		
	Sensor		CMOS
	Resolution	pixel	max. 2,592 x 1,944 (5 megapixels)
	Pilot laser		
	Wavelength	nm	635
	Output power	mW	1
	Lightning		
	LED lightning		
	Joystick		
	XYZ joystick with preset buttons		
	Sample		
	Sample weight including holder	kg	max. 5
	Holder for solid samples		yes
	Holder for liquid samples		optional
	Holder for gaseous samples		optional
<b>Safety</b>	Housing		laser class 1
	Laser protection window		adapted to applied LIBS laser wavelength
	Interlock		front door, detachable rear panel
<b>General</b>	Dimensions		
	Width x depth x height	mm	600 x 460 x 735
	Weight	kg	50
	Electric supply		
	Power supply	V DC	24
	Ports		
	Interlock		BNC
	Measurement signal output		SMA
	Gas input	mm	push-to-connect tube fitting, ∅ = 6
	Gas output	mm	push-to-connect tube fitting, ∅ = 6
	Exhaust flange	mm	tube, ∅ = 38
	Operating conditions		
	Temperature	°C	+ 10 ... + 40
	Relative humidity	%	max. 80

\* optional, other specifications possible

\*\* each axis

Subject to technical changes.